



Joseph P. Thompson & Associates, LLC
MUSCULOSKELETAL THERAPIES FOR ANIMALS
25 Kessel Court, Suite 103
MADISON, WI 53711 - USA
Tel : 608.274.7769
Fax : 262.364.2270
www.mtvet.com

Thermotex Therapeutic Pet Beds- Helping pets cope with pain.



Small Pet Bed

- ◆ Inside measures 11" x 13" x 4"
- ◆ Uses 110 volt, 35 watt power
- ◆ Polyester Blend interior / nylon exterior
- ◆ 3½" round foam cylinder enclosure
- ◆ For pets up to 12 lbs.

Thermotex pet beds are ideal for:

- ◆ post surgery recovery
- ◆ arthritis, bursitis, and strains
- ◆ muscles or joints that are causing pain.
- ◆ animals who need extra warmth.



Large Pet Bed

- ◆ Inside measures 12" x 20" x 4"
- ◆ Uses 110 volt, 35 watt power
- ◆ Polyester Blend interior / nylon exterior
- ◆ 3½" round foam cylinder enclosure
- ◆ For pets up to 24 lbs.

Thermotex pet beds have a Velcro opening for easy heat panel removal prior to washing the complete bed.

How does it work?

- * Thermotex Therapy Systems deliver pain relief by increasing the blood flow.
- * Increased blood flow is accomplished by using infrared energy absorption.
- * Infrared energy is absorbed by the body to a depth of 6 cm or 2.36 inches.
- * Infrared is a radiant energy source similar to the warmth felt by sunshine with no harmful rays.
- * As the infrared energy is absorbed, the blood flow and oxygen levels increase in the treated area.
- * This results in deep penetrating warmth.
- * This deep penetrating warmth delivers pain relief and assists with elimination of toxins.
- * Infrared is so safe, it is used to warm newborn babies in hospital.

Product Features

- * Inexpensive to operate: Thermotex Therapy Pet beds cost around 7 cents a day to operate for 24 hours.
- * Machine Washable Lining: Removable cloth liner is easily removed for washing.
- * No Burning: Thermotex Therapy Systems are thermostatically controlled and will never overheat.
- * Safety and simplicity of operation: Just plug the unit into any 110 volt outlet, set on low and leave it.
- * Uses 110 volt power; 35 watts.

Patents and Certifications

- * Patents pending
- * CSA / UL approved
- * ISO certified manufacturing standard 9001:2000