


**COMPANION ANIMAL
HYPERBARICS**

| | | Comments |
|------------------------------------|----------|---|
| Overall Length | 68.5" | |
| Internal Length | 58.5" | |
| Height @ Control Console | 58" | |
| Weight | 1182 lbs | By calculation |
| Chamber Diameter | 30" | |
| Chamber Greatest Width Dimension | 35.5" | At door closure, including hinge |
| Chamber Wall Thickness | 1/4" | Carbon steel |
| Internal Volume | 27.82 cf | Cubic Feet |
| Nominal Operating Pressure | 30 psi | |
| Maximum Allowable Working Pressure | 33 psi | Chamber PRV (Pressure Relief Valve) will exhaust at 35 psi |
| Hydrostatic Test Pressure | 48 psi | > 150 % of nominal pressure |
| Estimated Pressurization Time | < 10 min | To nominal working pressure |
| Estimated Depressurization Time | < 10 min | From nominal working pressure |
| | | All control components are readily available, "off the shelf" items |



CALL Joe Thompson – 866. 646. 8684